



Clever-Lite LS, Adult Walker with Seat and Push-Down Brakes



FIGURE A




FIGURE B

- Rear glide push-down brakes (Figure A)
- Sturdy 1" diameter anodized, extruded aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame
- Easily folds with patented dual-lever side-paddle release (Figure B)
- Soft, flexible backrest provides comfort and security while seated
- 5" fixed front wheels

PRODUCTS

PRODUCTS	Description
10233	Adult, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Width	18.5" inside hand grip 21" inside back legs 27" Overall
Depth	21" opened at base
Seat Dimensions	18.75" (W) x 11" (D)
Seat Height	19" - 25.5"
Handle Height	29.5" - 39"
Weight	11 lbs.
Carton Shipping Weight	15 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

PARTS

PARTS	Description
10243CPA	Seat Old Style (no frame), 1/ea
9500001232	Bushing, 1/ea
10325	Walker Tip Grey, 1/ea
10325B	Walker Tip Black, 1/ea
10109	Walker Wheel, 1/pr
12002PB	Push Button, 1/ea
10230SN	Seat w/Frame, 1/ea
10233RL	Rear Legs With Glide Brake, 1/pr
10230B	Basket, 1/ea
10230SF	Seat Frame, 1/ea