

Clinical Care Recliner



FIGURE A



FIGURE B



FIGURE D

- Ideal for use in dialysis, blood collection, respiratory and cardiac care
- Chair can be positioned in 4 positions—upright (Figure A), TV (Figure B), Full Recline (Figure C) and Trendelenburg
- Deep recline and Trendelenburg position can only be attained by an attendant
- Gas cylinder controls deep recline and Trendelenburg position
- Side panels "pop-off" for easy cleaning and maintenance
- Large, blow-molded side trays includes recess for cup
- 5" Casters (2 with locks) make it easy-to-move and secure recliner
- Gas cylinder release handle for quick back release
- Release handle can be operated from either side of the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Upholstered arm supports provide extra comfort in recline position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117 Standard

PRODUCTS	Description
D577-R	Rosewood, 1/cs
D577-J	Jade, 1/cs
D577-BR	Blue Ridge, 1/cs

SPECIFICATIONS	Description
Overall Dimensions	47.5" (H) x 30.5" (W) x 38" (D)
Width Between Arms	22"
Floor to Seat Height	26.5"
Seat Dimensions	21.5" (H) X 21.5" (W) X 20" (D)
Product Weight	115 lbs.
Carton Shipping Weight	135 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts

PARTS	Description
D577P-1003	1/ea
D577P-1002	1/ea
D577P-1001	1/ea
D577P-1000	1/ea
D577P-1005	1/ea
D577P-1060	1/ea
D574-1048	
D577BACK-BR	1/ea

PARTS	Description
D577BACK-J	1/ea
D577BACK-R	1/ea
D577SEAT-BR	1/ea
D577SEAT-J	1/ea
D577SEAT-R	1/ea
D577P-1010BR	1/ea
D577P-1010J	1/ea
D577P-1010R	1/ea
D574P-1035	1/bx
D577P-1012	1/ea
D577P-1013	1/ea
D577P-1014	1/ea
D577P-1075	1/bx
D577P-1026	1/bx
D574P-1025B	1/bx
D577P-1017	1/ea
D577P-1018BR	1/ea
D577P-1018J	1/ea
D577P -1018R	
D577P-1019	1/ea
D577P-1020	1/ea
D577P-1065	1/ea
D577P-1022	1/ea
D577P-1050	1/ea
D577P-1058	1/bx
D577P-1004BR	
D577P-1004J	
D577P-1004R	
D574P-1048	1/ea
D574P-1048-R	
D574P-1060	1/ea
D574P-1061	1/ea